


Gradient Plate Substitute

As you have already learned, random chance mutations occur in bacteria that enable them to resist antibiotics. In any given population of bacteria, there will be a few bacteria with naturally occurring mutations that provide some resistance to a particular antibiotic. We can select for those naturally occurring mutants using the gradient plate technique.


Two layers of agar are prepared: one layer is plain agar, the other layer contains streptomycin along with green food coloring so we can tell which layer is which. The plain agar is poured into a sterile petri plate and allowed to solidify at an angle. The streptomycin layer is poured on top. A concentration gradient of streptomycin will form. The darker green side has the highest concentration of streptomycin. The concentration decreases as you move across the plate to the low/no concentration side.

E. coli is spread all over the agar surface. After incubation, we see colonies growing in the high streptomycin concentration side of the plate.

Take Home Lesson: Why would you expect to find any streptomycin resistant mutants? **Random chance, spontaneous mutations.** What is the purpose of the streptomycin? **To select for resistant mutants.** Does streptomycin induce (cause) mutations? **No** What are some of the possible mechanisms a mutant might use to circumvent the streptomycin? **Review the lecture on antibiotics.**



Mutant colonies in high streptomycin concentration



Lots of *E. coli* growing in low/no streptomycin concentration