
LABORATORY BASICS FOR MICROBIOLOGY

How do we grow microorganisms?

Culture: the growth of bacteria (or any microorganisms) for the purpose of study

Culture medium: sterile mixture of nutrients (much like soup) for growing microorganisms




- 1. Liquid (broth)
- 2. Solid (agar)
 - a. Agar slants
 - b. Agar deeps
 - c. Agar plates



Questions:

- 1. What is the advantage of solid medium vs. liquid?
- 2. What is the use of an agar slant and agar deep?
- 3. Which container would provide the most surface area for growing bacteria?
- 4. What types of containers/culture vessels do we use?
- 5. In which container would you put liquid medium?
- 6. In which container would you put solid medium?



Write one use for each of the followings:

Test tube with soft agar:

Test tubes with screw cap

Test tube with very little agar at the bottom:

Test tubes with colored broth:

Test tubes with different color tops:

Test tubes with durham tube:

Empty test tube with no cap:

- 7. Which container would you GUESS to be the best for bacteria which cannot tolerate air?
- 8. What would be the DISADVANTAGE of the screw cap?

Microbiological test tubes with various closures:

- A. Metal cap (notice tiny holes) pulls off:
- B. Screw cap
- C. Plastic cap (no holes/pulls off):
- D. Open tube (No cap)

How do we transfer liquids and microorganisms from one container to another

Transfer liquids:

Pipettes (glass/plastic)

Pipetters/Pipette aid

Transfer microorganisms:

Spreader

Inoculating loop

Inoculating needle

Spreader

Sterile swabs

Questions:

Queblions.	
A: 10.0ml B: 1.0ml	9. What is the <i>total</i> volume you can transfer with pipette A ?
DNT at the P	10. What is the <i>total</i> volume you can transfer with pipette B ?
	11. How does the bracketed section \mathbf{X} of pipette \mathbf{A} differ from that of pipette \mathbf{B} ?
	12. What is the difference between each line marked by section Y of pipette A vs. pipette B?
	13. Which of the two pipettes, A or B should be used to accurately measure 0.1 ml?
5 Y	

