

Microbiology Lab Experiment Changes

Experiment #: 3-7

Title: Acid-Fast Stain

Live Organisms: *Mycobacterium phlei*, *E. coli*

Changes:

Modified procedure:

1. Make *E. coli* smear slightly cloudy as usual. For the *M. phlei* smear, just start mixing right away. Break up clumps as best as you can. Smear will be clumpy.
2. Put **heat fixed** smear onto "staining hot plate" in fume hood.
3. Place a filter paper over smear. Saturate filter paper with carbol fuchsin stain.
4. Allow stain to sit approximately 5 minutes.
5. Get a paper towel to carry slide to sink. Throw filter paper into trash. Rinse slide well.
6. Decolorize with Acid-Alcohol: 1-2 drops for *M. phlei* and five+ drops for *E. coli*.
7. Counterstain with Methylene Blue for 1-2 minutes.
8. Rinse with water. Blot dry. You may wipe **BOTTOM** of slide with lens cleaning fluid and tissue.

Each student will make one slide of each organism.

Take Home Lesson: You need to know the reagents and their function in the staining procedure. Why do we use the acid-fast stain? On what organisms is it used? Be able to recognize that an organism is acid-fast or non-acid-fast based on its appearance under the microscope.