Microbiology Lab Experiment Changes

Experiment #: 4-1, 4-2, 4-4, 4-6

Title: Selective and Differential Media

Live Organisms: Enterobacter aerogenes, Escherichia coli, Bacillus subtilis,

Staphylococcus aureus, Staphylococcus epidermidis, Providencia rettgeri

Changes: Procedure:

Work in groups (2-3 students). Each person in your group should do at least

one plate.

1. Inoculate 4 different bacteria onto each kind of medium.

MSA	PEA	MAC	EMB
 S. aureus S. epidermidis B. subtilis E. coli or aerogenes or rettgeri 	 S. aureus S. epidermidis B. subtilis E. coli or aerogenes or rettgeri 	 E. coli E. aerogenes P. rettgeri S. aureus or S. epidermidis or B. subtilis 	 E. coli E. aerogenes P. rettgeri S. aureus or epidermidis or subtilis

- 2. The above combinations demonstrate the selective and / or differential nature of each plate.
- 3. For the sake of conformity and reducing color confusion, we have adopted the following color conventions for bacterial growth on the following agar plates:

On MacConkey agar: lactose fermenters (LF) = "pink/magenta" non-lactose fermenters (NLF) = "colorless" (i.e., greyish, whiteish) On EMB agar: LF = "dark purple" NLF = "purplish-blue"

4. Next lab period record and draw results. Note color of medium, color of bacterial growth, presence of precipitate, and the absence or inhibition of growth.

Take Home Lesson: Know the principles described in the lab manual for each type of medium. Distinguish between selective and differential media. Understand that a particular medium may be both selective and differential. Describe the selective and differential nature of: EMB, MAC, PEA and MSA. What are the components in each medium and what is each component's function with respect to the medium's selective and / or differential properties?

Print and study the Selective and Differential Media handout from the web site.